

# Specimen Collection and Submission Instructions

## Mumps IgM, EIA (Version 1)

Page 1 of 2

Specimen Type	Collection Time	Collection Frequency	Collection Procedures	Transport Media	Shipping & Handling (S&H)**
Blood	As soon as possible upon suspicion of disease.	NA	<ul style="list-style-type: none"> <li>Collect at least 4 mL of blood in a red top or red-gray top tube.</li> <li><b>Minimum serum volume:</b> 1 mL.</li> </ul>	NA	<ul style="list-style-type: none"> <li><b>Transport device:</b> Red top or red-gray top tube.</li> <li><b>Transport:</b> Refrigerate and ship cold (2-8°C) on ice packs to arrive at WAPHL during business hours within one week of collection. If previously frozen, ship on dry ice. Ship as Category B.</li> <li><b>Storage:</b> Refrigerate at 2-8°C up to one week after collection. For long-term storage, freeze at -20°C.</li> <li><b>Rejection Criteria:</b> <ul style="list-style-type: none"> <li>Specimen collected too early (before rash onset).</li> <li>Missing requisition form.</li> <li>No submitter information on form.</li> <li>Specimen is excessively hemolyzed or lipemic.</li> <li>Missing collection date.</li> <li>Missing rash onset date.</li> <li>Missing immunization status of patient.</li> <li>Please see general rejection criteria on page 2.</li> </ul> </li> </ul>
Serum*					

## Specimen Collection and Submission Instructions

### Mumps IgM, EIA (Version 1)

Page 2 of 2

#### **General Rejection Criteria (for additional details, see S&H)**

- Unaccepted specimen type.
- Not refrigerated or frozen properly.
- Insufficient specimen volume.
- Leaky specimen.
- Did not obtain LHJ approval.
- Failure to follow specific S&H requirements.

#### **NOTE:**

Specimens will not be processed until ALL of the following information are provided on requisition form:

- Patient name, patient's county of residence, second identifier, and date of onset.
- Two matching identifiers on both the tube and the laboratory requisition form.
- Specimen type, date of collection, and test requested.
- Submitter name, address, and telephone/FAX numbers.

#### **Collection kit availability**

Not available.

#### **Comments**