

# Specimen Collection and Submission Instructions

## Antimicrobial Susceptibility Testing (Version 1)

Page 1 of 2

Specimen Type	Collection Time	Collection Frequency	Collection Procedures	Transport Media	Shipping & Handling (S&H)**
Isolate	NA	NA	NA	Appropriate growth media	<ul style="list-style-type: none"> <li>• <b>Transport device:</b> Screw-cap tube or agar plate (seal with parafilm, if available).</li> <li>• <b>Rejection Criteria:</b> <ul style="list-style-type: none"> <li>○ Mixed culture.</li> <li>○ Fungal contamination.</li> <li>○ Expired media.</li> <li>○ Broken transport device.</li> <li>○ Leaking specimen.</li> </ul> </li> <li>• <b>Transport:</b> Ship ASAP at ambient temperature. Ship as Category B.</li> <li>• <b>Storage:</b> Store at ambient temperature.</li> </ul>

# Specimen Collection and Submission Instructions

## Antimicrobial Susceptibility Testing (Version 1)

Page 2 of 2

### **General Rejection Criteria (for additional details, see S&H)**

- Failure to follow specific S&H requirements.

**NOTE:**

Specimens will not be processed until ALL of the following information are provided on requisition form:

- Patient name, patient's county of residence, and second identifier.
- Two matching identifiers on both the tube and the laboratory requisition form.
- Specimen type, date of collection, and test requested.
- Submitter name, address, and telephone/FAX numbers.
- Reason for testing at WAPHL.

### **Collection kit availability**

Not available

### **Comments**