

Facility/MTS # _____ Date _____

Y N WAC 246-338-020 LICENSURE

— — The Medical Test Site has a current license appropriate for the services provided.

Y N NA WAC 246-338-050 PROFICIENCY TESTING

— — — (1) All licensed medical test sites, excluding those granted a certificate of waiver, must:

— — — (a) Comply with federal proficiency testing requirements listed in 42 CFR Part 493-Laboratory Requirements, Subparts H and I;

Y N NA WAC 246-338-060 PERSONNEL

— — — (1) Medical test site owners must:

— — — (a) Have a director responsible for the overall technical supervision and management of the test site personnel including oversight of the performance of test procedures and reporting of test results;

— — — (b) Have technical personnel, competent to perform tests and report test results; and

— — — (c) Meet the standards for personnel qualifications and responsibilities in compliance with federal regulation, as listed in 42 CFR Part 493 Subpart M-Personnel for Moderate and High Complexity Testing.

— — — (3) Medical test site directors must:

— — — (a) Establish and approve policies for:

— — — (i) Performing, recording, and reporting of tests;

— — — (ii) Maintaining an ongoing quality assurance program;

— — — (iii) Supervision of testing; and

— — — (iv) Compliance with chapter 70.42 RCW and this chapter;

— — — (b) Evaluate, verify, and document the following related to technical personnel:

— — — (i) Education, experience, and training in test performance and reporting test results;

— — — (ii) Sufficient numbers to cover the scope and complexity of the services provided;

— — — (iii) Access to training appropriate for the type and complexity of the test site services offered; and

— — — (iv) Maintenance of competency to perform test procedures and report test results;

— — — (c) Be present, on call, or delegate the duties of the director to an on-site technical person during testing.

Y N NA WAC 246-338-070 RECORDS

— — — Medical test sites must maintain records as described in this section.

— — — (1) REQUISITIONS must include the following information, in written or electronic form:

— — — (a) Patient name, identification number, or other method of specimen identification;

— — — (b) Name or other suitable identifier of the authorized person ordering the test;

— — — (c) Date of specimen collection, and time, if appropriate;

— — — (d) Source of specimen, if appropriate;

— — — (e) Type of test ordered;

— — — (f) Sex and age of the patient, if appropriate; and

Facility/MTS # _____ Date _____

Y N NA

WAC 246-338-070 RECORDS

- (g) For cytology and histopathology specimens:
 (i) Pertinent clinical information; and
 (ii) For Pap smears:
 (A) Date of last menstrual period; and
 (B) Indication whether the patient has history of cervical cancer or its precursors.
- (2) TEST RECORD SYSTEMS must:
 (a) Consist of instrument printouts, worksheets, accession logs, corrective action logs, and other records that ensure reliable identification of patient specimens as they are processed and tested to assure that accurate test results are reported; and
 (b) Include:
 (i) The patient's name or other method of specimen identification;
 (ii) The date the specimen was received, and time, if appropriate;
 (iii) The reason for specimen rejection or limitation;
 (iv) The date of specimen testing; and
 (v) The identification of the personnel who performed the test.
- (3) TEST REPORTS must:
 (a) Be maintained in a manner permitting identification and reasonable accessibility;
 (b) Be released only to authorized persons or designees;
 (c) Include the name and address of the medical test site, or where applicable, the name and address of each medical test site performing each test;
 (d) Include:
 (i) Date reported;
 (ii) Time reported, if appropriate;
 (iii) Any information regarding specimen rejection or imitation; and
 (iv) Name of the test performed, test result, and units of measurement, if applicable.
- (4) CYTOLOGY REPORTS must:
 (a) Distinguish between unsatisfactory specimens and negative results;
 (b) Provide narrative descriptions for any abnormal results, such as the Bethesda system of terminology as published in the Journal of the American Medical Association, 1989, Volume 262, pages 931-934; and
 (c) Include the signature or initials of the technical supervisor, or an electronic signature authorized by the technical supervisor, for nongynecological preparations and gynecological preparations interpreted to be showing reactive or reparative changes, atypical squamous or glandular cells of undetermined significance, or to be in the premalignant (dysplasia, cervical intraepithelial neoplasia or all squamous intraepithelial neoplasia lesions including human papillomavirus-associated changes) or malignant category.
- (5) HISTOPATHOLOGY REPORTS must include the signature or initials of the technical supervisor or an electronic signature authorized by the technical supervisor on all reports.
- (6) CYTOGENETICS REPORTS must:
 (a) Use appropriate nomenclature on final reports;
 (b) Include the number of cells counted and karyotyped; and
 (c) Include an interpretation of the karyotypes findings.
- (7) If a specimen is referred to another laboratory for testing, the medical test site must:

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-070 RECORDS

- — — (a) Report the essential elements of the referred test results without alterations that could affect the clinical interpretation of the results; and
- — — (b) Retain or be able to produce an exact duplicate of each testing report from the referral laboratory.
- — — (8) The medical test site must retain records, slides, and tissues as described in Table 070-1.

Table 070-1 Record/Slide/Tissue Retention Schedule

	Two Years	Five Years	Ten Years
(a) General Requirements for all Laboratory Specialties	<ul style="list-style-type: none"> ▪ Test requisitions or equivalent; ▪ Test records; ▪ Test reports; ▪ Quality control records; ▪ Quality assurance records; ▪ Proficiency testing records; ▪ Hard copy of report, or ability to reproduce a copy, for all specimens referred for testing; and ▪ Discontinued procedures for all specialty areas 		
(b) Transfusion Services*		<ul style="list-style-type: none"> ▪ Test requisitions or equivalent; ▪ Test records; ▪ Test reports; ▪ Quality control records; and ▪ Quality assurance records 	
(c) Cytology		<ul style="list-style-type: none"> ▪ All cytology slides, from date of examination of the slide 	<ul style="list-style-type: none"> ▪ All cytology reports
(d) Histopathology	<ul style="list-style-type: none"> ▪ Specimen blocks, from date of examination 		<ul style="list-style-type: none"> ▪ All histopathology reports; and ▪ Stained slides, from date of examination of the slide
(e) Histopathology – Tissues	<ul style="list-style-type: none"> ▪ Retain remnants of tissue specimens in an appropriate preserved state until the portions submitted for microscopic examination have been examined and diagnosed 		
(f) Instrument/method Validation Studies	<ul style="list-style-type: none"> ▪ For life of instrument/method plus two years 		

*Must be retained for no less than five years in accordance with 21 CFR Part 606, Subpart I

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-080 QUALITY ASSURANCE

- | | | | |
|---|---|---|---|
| — | — | — | Each medical test site performing moderate complexity (including PPMP) or high complexity testing, or any combination of these tests, must establish and follow written policies and procedures for a comprehensive quality assurance program. The quality assurance program must be designed to monitor and evaluate the ongoing and overall quality of the total testing process (preanalytic, analytic, postanalytic). The medical test site's quality assurance program must evaluate the effectiveness of its policies and procedures; identify and correct problems; assure the accurate, reliable, and prompt reporting of test results; and assure the adequacy and competency of the staff. As necessary, the medical test site must revise policies and procedures based upon the results of those evaluations. The medical test site must meet the standards as they apply to the services offered, complexity of testing performed and test results reported, and the unique practices of each testing entity. All quality assurance activities must be documented. |
| | | | (1) The medical test site must establish and implement a written quality assurance plan, including policies and procedures, designed to: |
| — | — | — | (a) Monitor, evaluate, and review quality control data, proficiency testing results, and test results, including biannual verification of: |
| — | — | — | (i) Accuracy of test results for tests that are not covered by proficiency testing; and |
| — | — | — | (ii) Relationship between test results when the medical test site performs the same test on different instruments or at different locations within the medical test site; |
| — | — | — | (b) Identify and correct problems; |
| — | — | — | (c) Establish and maintain accurate, reliable, and prompt reporting of test results; |
| — | — | — | (d) Verify all tests performed and reported by the medical test site conform to specified performance criteria in quality control under WAC 246-338-090; and |
| — | — | — | (e) Establish and maintain the adequacy and competency of the technical personnel. |
| | | | (2) The quality assurance plan must include mechanisms or systems to: |
| — | — | — | (a) Establish and apply criteria for specimen acceptance and rejection; |
| — | — | — | (b) Notify the appropriate individuals as soon as possible when test results indicate potential life-threatening conditions; |
| — | — | — | (c) Assess problems identified during quality assurance reviews and discuss them with the appropriate staff; |
| — | — | — | (d) Evaluate all test reporting systems to verify accurate and reliable reporting, transmittal, storage, and retrieval of data; |
| — | — | — | (e) Document all action taken to identify and correct problems or potential problems; |
| — | — | — | (f) Issue corrected reports when indicated; |
| — | — | — | (g) Provide appropriate instructions for specimen collection, handling, preservation, and transportation; and |
| — | — | — | (h) Provide clients updates of testing changes that would affect test results or the interpretation of test results. |
| — | — | — | (3) The medical test site must establish criteria for and maintain appropriate documentation of any remedial action taken in response to quality control, quality assurance, personnel, proficiency testing, and transfusion reaction investigations. |

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-080 QUALITY ASSURANCE

- (4) The medical test site must have a system in place to assure:
- — — (a) All complaints and problems reported to the medical test site are documented and investigated when appropriate; and
- — — (b) Corrective actions are instituted as necessary.
- (5) The owner must:
- — — (a) Maintain adequate space, facilities, and essential utilities for the performance and reporting of tests;
- — — (b) Establish, post, and observe safety precautions to ensure protection from physical, chemical, biochemical, and electrical hazards and biohazards; and
- — — (c) Establish and implement policies and procedures for infectious and hazardous medical wastes consistent with local, state, and federal authorities.
- (6) Information that must be available to authorized persons ordering or utilizing the test results includes:
- — — (a) A list of test methods, including performance specifications;
- — — (b) Reference ranges; and
- — — (c) Test method limitations
- (7) If the medical test site refers specimens to another site for testing, the site to which specimens are referred must have a valid medical test site license or meet equivalent requirements as determined by HCFA.

Y N NA WAC 246-338-090 QUALITY CONTROL

- The medical test site must use quality control procedures, providing and assuring accurate and reliable test results and reports, meeting the requirements of this chapter.
- (1) The medical test site must have written procedures and policies available in the work area for:
- — — (a) Analytical methods used by the technical personnel including:
- (i) Principle;
- (ii) Specimen collection and processing procedures;
- (iii) Equipment/reagent/supplies required;
- (iv) Preparation of solutions, reagents, and stains;
- (v) Test methodology;
- (vi) Quality control procedures;
- (vii) Procedures for reporting results (normal, abnormal, and critical values);
- (viii) Reference range;
- (ix) Troubleshooting guidelines - limitations of methodology;
- (x) Calibration procedures; and
- (xi) Pertinent literature references; and
- — — (b) Alternative or backup methods for performing tests including the use of a reference facility if applicable.
- (2) The medical test site must establish written criteria for and maintain appropriate documentation of:
- — — (a) Temperature-controlled spaces and equipment;
- — — (b) Preventive maintenance activities;
- — — (c) Equipment function checks;
- — — (d) Procedure calibrations; and
- — — (e) Method/instrument validation procedures.
- (3) The medical test site must maintain documentation of:
- — — (a) Expiration date, lot numbers, and other pertinent information for:
- (i) Reagents;

Facility/MTS # _____ Date _____

Y N NA

WAC 246-338-090 QUALITY CONTROL

- | | | | |
|---|---|---|---|
| — | — | — | (ii) Solutions; |
| — | — | — | (iii) Culture media; |
| — | — | — | (iv) Controls; |
| — | — | — | (v) Calibrators; |
| — | — | — | (vi) Standards; |
| — | — | — | (vii) Reference materials; and |
| — | — | — | (viii) Other testing materials; and |
| — | — | — | (b) Testing of quality control samples. |
| | | | (4) For quantitative tests , the medical test site must perform quality control as follows: |
| — | — | — | (a) Include two reference materials of different concentrations each day of testing unknown samples, if these reference materials are available; or |
| — | — | — | (b) Have an equivalent mechanism to assure the quality, accuracy, and precision of the test if reference materials are not available. |
| | | | (5) For qualitative tests , the medical test site must perform quality control as follows: |
| — | — | — | (a) Use positive and negative reference material each day of testing unknown samples; or |
| — | — | — | (b) Have an equivalent mechanism to assure the quality, accuracy, and precision of the test if reference materials are not available. |
| | | | (6) The medical test site must: |
| — | — | — | (a) Use materials within their documented expiration date; |
| — | — | — | (b) Not interchange components of kits with different lot numbers, unless specified by the manufacturer; |
| — | — | — | (c) Determine the statistical limits for each lot number of unassayed reference materials through repeated testing; |
| | | | (d) Use the manufacturer's reference material limits for assayed material, provided they are: |
| — | — | — | (i) Verified by the medical test site; and |
| — | — | — | (ii) Appropriate for the methods and instrument used by the medical test site; |
| — | — | — | (e) Make reference material limits readily available; |
| — | — | — | (f) Report patient results only when reference materials are within acceptable limits; and |

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 QUALITY CONTROL

— — — (g) Comply with general quality control requirements as described in Table 090-1, unless otherwise specified in subsection (9)(a) through (l) of this section.

Table 090-1 General Quality Control Requirements

	Control Material	Frequency
(a) Each batch or shipment of reagents, discs, antisera, and identification systems	• Appropriate control materials for positive and negative reactivity	• When prepared or opened, unless otherwise specified
(b) Each batch or shipment of stains	• Appropriate control materials for positive and negative reactivity	• When prepared or opened; and • Each day of use, unless otherwise specified
(c) Fluorescent stains	• Appropriate control materials for positive and negative reactivity	• Each time of use, unless otherwise specified
(d) Quality control for each specialty and subspecialty	• Appropriate control materials; or • Equivalent mechanism to assure the quality, accuracy, and precision of the test if reference materials are not available	• At least as frequently as specified in this section; • More frequently if recommended by the manufacturer of the instrument or test procedure; or • More frequently if specified by the medical test site
(e) Direct antigen detection systems without procedural controls	• Positive and negative controls that evaluate both the extraction and reaction phase	• Each batch, shipment, and new lot number; and • Each day of use

Y N NA WAC 246-338-090 QUALITY CONTROL

— — — (7) The medical test site must perform, when applicable:
(a) Calibration and calibration checks for **moderate complexity testing** as described in Table 090-2;

Table 090-2 Calibration and Calibration Checks - Moderate Complexity Testing

	Calibration Material	Frequency
CALIBRATION	• Calibration material appropriate for methodology according to manufacturer's instructions	• Initial on-site installation/implementation of instrument/method; • At the frequency recommended by the manufacturer; • When controls show trends, shifts, or are out of limits and other corrective action has not fixed the problem.
CHECK CALIBRATION	• Assayed material appropriate for methodology	• At least every six months.

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 QUALITY CONTROL

— — — (b) Calibration and calibration verification for **high complexity testing** as described in Table 090-3;

Table 090-3 Calibration and Calibration Checks - High Complexity Testing

	Calibration Material	Frequency
CALIBRATION	<ul style="list-style-type: none"> Calibration materials appropriate for methodology 	<ul style="list-style-type: none"> Initial on-site installation/implementation of instrument/method; At the frequency recommended by the manufacturer; and Whenever calibration verification fails to meet the medical test site's acceptable limits for calibration verification.
CALIBRATION VERIFICATION	<ul style="list-style-type: none"> Use assayed material, if available, at the lower, mid-point, and upper limits of procedure's reportable range; or Demonstrate alternate method of assuring accuracy at the lower, mid-point, and upper limits of procedure's reportable range 	<ul style="list-style-type: none"> At least every six months; When there is a complete change of reagents (<i>i.e.</i>, new lot number or different manufacturer) is introduced; When major preventive maintenance is performed or there is a replacement of critical parts of equipment; or When controls are outside of the medical test site's acceptable limits or exhibit trends.

Y N NA WAC 246-338-090 QUALITY CONTROL

— — — (c) Validation for **moderate complexity testing** by verifying the following performance characteristics when the medical test site introduces a new procedure classified as moderate complexity:

- (i) Accuracy;
- (ii) Precision; and
- (iii) Reportable range of patient test results;

(d) Validation for **high complexity testing**:

- (i) When the medical test site introduces a new procedure classified as high complexity;
- (ii) For each method that is developed in-house, is a modification of the manufacturer's test procedure, or is an instrument, kit or test system that has not been cleared by FDA; and
- (iii) By verifying the following performance characteristics:
 - (A) Accuracy;
 - (B) Precision;
 - (C) Analytical sensitivity;
 - (D) Analytical specificity to include interfering substances;
 - (E) Reference ranges (normal values);
 - (F) Reportable range of patient test results; and
 - (G) Any other performance characteristic required for test performance.

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 QUALITY CONTROL

(8) When patient values are above the maximum or below the minimum calibration point or the reportable range, the medical test site must:

—	—	—	(a) Report the patient results as greater than the upper limit or less than the lower limit or an equivalent designation; or
---	---	---	--

—	—	—	(b) Use an appropriate procedure to rerun the sample allowing results to fall within the established linear range.
---	---	---	--

—	—	—	(9) The medical test site must perform quality control procedures as described for each specialty and subspecialty in (a) through (l) of this subsection.
---	---	---	---

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL— — — (a) **Chemistry:** Perform quality control procedures for chemistry as described in Table 090-4.**Table 090-4 Quality Control Procedures - Chemistry**

Subspecialty/Test		Qualitative		Quantitative	
	Control Material		Frequency	Control Material	Frequency
Routine Chemistry	• Positive and negative reference material	•	Each day of use	• Two levels of reference material in different concentrations	• Each day of use
Toxicology					
• GC/MS for drug screening	• Analyte-specific control	•	With each run of patient specimens	• Analyte-specific control	• With each analytical run
• Urine drug screen	• Positive control containing at least one drug representative of each drug class to be reported; must go through each phase of use including extraction	•	With each run of patient specimens		
Urinalysis					
• Non-waived instrument				• Two levels of control material	• Each day of use
• Refractometer for specific gravity				• Calibrate to zero with distilled water • One level of control material	• Each day of use
Blood Gas Analysis				• Two-point calibration and one reference material	• Each 8 hours of testing
				• One-point calibration or one reference material, or	• Each time patient sample is tested, unless automated instrument internally verifies calibration every 30 minutes
				• Another calibration and reference material schedule, approved by the department	
Electrophoresis	• One control containing fractions representative of those routinely reported in patient specimens	•	In each electrophoretic cell	• One control containing fractions representative of those routinely reported in patient specimens	• In each electrophoretic cell

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL**(b) Hematology:**

- | | | | |
|---|---|---|--|
| — | — | — | (i) Run patient and quality control samples in duplicate for manual cell counts; |
| — | — | — | (ii) If reference material is unavailable, document the mechanism used to assure the quality, accuracy, and precision of the test; and |
| — | — | — | (iii) Perform quality control procedures for hematology as described in Table 090-5. |

Table 090-5 Quality Control Procedures - Hematology

	Control Material	Frequency
Automated	• Two levels of reference material in different concentrations	• Every 8 hours that patient samples are tested
Manual Blood Counts	• One level of reference material	• Every 8 hours that patient samples are tested
Qualitative Tests	• Positive and negative reference material	• Each day of testing

Y N NA WAC 246-338-090 (9) QUALITY CONTROL**(c) Coagulation:**

- | | | | |
|---|---|---|---|
| — | — | — | (i) Run patient and quality control samples in duplicate for manual coagulation test (tilt tube); |
| — | — | — | (ii) If reference material is unavailable, document the mechanism used to assure the quality, accuracy and precision of the test; and |
| — | — | — | (iii) Perform quality control procedures for coagulation as described in Table 090-6. |

Table 090-6 Quality Control Procedures - Coagulation

	Control Material	Frequency
Automated	• Two levels of reference material in different concentrations	• Every 8 hours that patient samples are tested; and • Each time reagents are changed
Manual Tilt Tube Method	• Two levels of reference material in different concentrations	• Every 8 hours that patient samples are tested; and • Each time reagents are changed

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL**(d) General Immunology:**

- | | | | |
|---|---|---|--|
| — | — | — | (i) Employ reference materials for all test components to ensure reactivity; |
| — | — | — | (ii) Report test results only when the predetermined reactivity pattern of the reference material is observed; and |
| — | — | — | (iii) Perform quality control procedures for general immunology as described in Table 090-7. |

Table 090-7 Quality Control Procedures - General Immunology

	Control Material	Frequency
Serologic tests on unknown specimens	• Positive and negative reference material	• Each day of testing
Moderate complexity kits with procedural (internal) controls	• Positive and negative reference material (external controls) • Procedural (internal) controls	• When kit is opened • Each time patient sample is tested

Y N NA WAC 246-338-090 (9) QUALITY CONTROL**(e) Syphilis Serology:**

- | | | | |
|---|---|---|---|
| — | — | — | (i) Use equipment, glassware, reagents, controls, and techniques that conform to manufacturer's specifications; |
| — | — | — | (ii) Employ reference materials for all test components to ensure reactivity; and |
| — | — | — | (iii) Perform serologic tests on unknown specimens concurrently with a positive serum reference material with known titer or graded reactivity and a negative reference material. |

(f) Microbiology:

- | | | | |
|---|---|---|--|
| — | — | — | (i) Have available and use: |
| — | — | — | (A) Appropriate stock organisms for quality control purposes; and |
| — | — | — | (B) A collection of slides, photographs, gross specimens, or text books for reference sources to aid in identification of microorganisms; |
| — | — | — | (ii) Document all steps (reactions) used in the identification of microorganisms on patient specimens; |
| — | — | — | (iii) For antimicrobial susceptibility testing: |
| — | — | — | (A) Record zone sizes or minimum inhibitory concentration for reference organisms; and |
| — | — | — | (B) Zone sizes or minimum inhibitory concentration for reference organisms must be within established limits before reporting patient results; and |
| — | — | — | (C) Perform quality control on antimicrobial susceptibility testing media as described in Table 090-9; |
| — | — | — | (iv) For noncommercial media, check each batch or shipment for sterility, ability to support growth and, if appropriate, selectivity, inhibition, or biochemical response; |

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

(v) For commercial media:

- — — (A) Verify that the product insert specifies that the quality control checks meet the requirements for media quality control as outlined by the National Committee for Clinical Laboratory Standards (NCCLS), Quality Assurance for Commercially Prepared Microbiological Culture Media-Second Edition; Approved Standard (1996);
- — — (B) Keep records of the manufacturer's quality control results;
- — — (C) Document visual inspection of the media for proper filling of the plate, temperature or shipment damage, and contamination before use; and
- — — (D) Follow the manufacturer's specifications for using the media; and

(vi) For microbiology subspecialties:

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

- — — (A) **Bacteriology:** Perform quality control procedures for bacteriology as described in Tables 090-8 and 090-9.

Table 090-8 Quality Control Procedures - Bacteriology

	Control Material	Frequency
Reagents, disks, and identification systems	<ul style="list-style-type: none"> Positive and negative reference organisms, unless otherwise specified 	<ul style="list-style-type: none"> Each batch, shipment and new lot number unless otherwise specified
Stains, unless otherwise specified; DNA probes; catalase; coagulase; beta-lactamase; and oxidase reagents	<ul style="list-style-type: none"> Positive and negative reference organisms 	<ul style="list-style-type: none"> Each batch, shipment and new lot number; and Each day of use
Fluorescent stains	<ul style="list-style-type: none"> Positive and negative reference organisms 	<ul style="list-style-type: none"> Each batch, shipment and new lot number; and Each time of use
Gram and acid-fast stains, bacitracin, optochin, ONPG, X and V disks or strips	<ul style="list-style-type: none"> Positive and negative reference organisms 	<ul style="list-style-type: none"> Each batch, shipment and new lot number; and Each week of use
Direct antigen detection systems without procedural controls	<ul style="list-style-type: none"> Positive and negative controls that evaluate both the extraction and reaction phase 	<ul style="list-style-type: none"> Each batch, shipment and new lot number; and Each day of use
Moderate complexity test kits with procedural (internal) controls	<ul style="list-style-type: none"> Positive and negative reference material (external) controls Procedural (internal) controls 	<ul style="list-style-type: none"> Each batch, shipment and new lot number Each time patient sample is tested
Antisera	<ul style="list-style-type: none"> Positive and negative reference material 	<ul style="list-style-type: none"> Each batch, shipment and new lot number; and Each month of use

Facility/MTS # _____ Date _____

Table 090-9 Quality Control Procedures - Bacteriology – Media for Antimicrobial Susceptibility Testing

	Control Material	Frequency
Check each new batch of media and each new lot of antimicrobial disks or other testing systems (MIC)	<ul style="list-style-type: none"> Approved reference organisms (ATCC organisms) 	<ul style="list-style-type: none"> Before initial use and each day of testing; or May be done weekly if the medical test site can meet the quality control requirements for antimicrobial disk susceptibility testing as outlined by NCCLS Performance Standards for Antimicrobial Disk Susceptibility Tests-Seventh Edition; Approved Standard (2000)

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

___ ___ ___ **(B) Mycobacteriology:** Perform quality control procedures for mycobacteriology as described in Table 090-10.

Table 090-10 Quality Control Procedures - Mycobacteriology

	Control Material	Frequency
Iron uptake test	<ul style="list-style-type: none"> Acid-fast organism that produces a positive reaction and with an organism that produces a negative reaction 	<ul style="list-style-type: none"> Each day of use
All other reagents or test procedures used for mycobacteria identification unless otherwise specified	<ul style="list-style-type: none"> Acid-fast organism that produces a positive reaction 	<ul style="list-style-type: none"> Each day of use
DNA probes	<ul style="list-style-type: none"> Organisms that produce positive and negative reactions 	<ul style="list-style-type: none"> Each day of use
Acid-fast stains	<ul style="list-style-type: none"> Acid-fast organism that produces a positive reaction 	<ul style="list-style-type: none"> Each week of use
Fluorochrome acid-fast stains	<ul style="list-style-type: none"> Organisms that produce positive and negative reactivity 	<ul style="list-style-type: none"> Each week of use
Susceptibility tests performed on <i>Mycobacterium tuberculosis</i> isolates	<ul style="list-style-type: none"> Strain of <i>M. tb</i> susceptible to all antimycobacterial agents used 	<ul style="list-style-type: none"> Each week of use

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

— — — (C) **Mycology:** Perform quality control procedures for mycology as described in Table 090-11.

Table 090-11 Quality Control Procedures - Mycology

	Control Material	Frequency
Auxanographic medium for nitrate assimilation: nitrate reagent	• Peptone control	• Each day of use
Susceptibility tests: each drug NOTE: Establish control limits and criteria for acceptable control results prior to reporting patient results	• One control strain that is susceptible to the drug	• Each day of use
Acid-fast stains	• Organisms that produce positive and negative reactions	• Each week of use
Reagents for biochemical and other identification test procedures	• Organism that produces a positive reaction	• Each week of use
Commercial identification systems utilizing 2 or more substrates	• Organisms that verify positive and negative reactivity of each media type	• Each batch or shipment and each lot number

Y N NA WAC 246-338-090 (9) QUALITY CONTROL**(D) Parasitology:**

(I) Have available and use:

- — — • Reference collection of slides or photographs and, if available, gross specimens for parasite identification; and
- — — • Calibrated ocular micrometer for determining the size of ova and parasites, if size is a critical parameter.

(II) Check permanent stains each month of use with reference materials.

E) Virology:

(I) Have available:

- — — • Host systems for isolation of viruses; and
- — — • Test methods for identification of viruses that cover the entire range of viruses that are etiologically related to the clinical diseases for which services are offered; and
- — — (II) Simultaneously culture uninoculated cells or cell substrate as a negative control when performing virus identification.

(g) **Histopathology:** Include a control slide of known reactivity with each slide or group of slides for differential or special stains and document reactions.**(h) Cytology:**

(i) Processing Specimens:

- — — (A) Stain all gynecological smears using a Papanicolaou or a modified Papanicolaou staining method;

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

- | | | | |
|---|---|---|--|
| — | — | — | (B) Have methods to prevent cross-contamination between gynecologic and nongynecologic specimens during the staining process; and |
| — | — | — | (C) Stain nongynecological specimens that have a high potential for cross-contamination separately from other nongynecological specimens, and filter or change the stains following staining. |
| | | | (ii) Performing Specimen Examinations: |
| — | — | — | (A) All cytology preparations must be evaluated on the premises of the medical test site; |
| — | — | — | (B) Technical personnel must examine, unless federal law and regulation specify otherwise, no more than one hundred cytological slides by nonautomated microscopic technique in a twenty-four-hour period and in no less than an eight-hour work period; |
| — | — | — | (C) Previously examined negative, reactive, reparative, atypical, premalignant or malignant gynecological cases and previously examined nongynecologic cytology preparations and tissue pathology slides examined by a technical supervisor are not included in the one hundred slide limit; |
| — | — | — | (D) Each slide preparation technique (automated, semi-automated or liquid based) which results in cell dispersion over one-half or less of the total available slide area and which is examined by nonautomated microscopic technique must be counted as one-half slide; and |
| — | — | — | (E) Records of the total number of slides examined by each individual at all sites during each twenty-four-hour period must be maintained. |
| | | | (iii) Establish and implement a quality assurance program that ensures: |
| — | — | — | (A) There is criteria for submission of material; |
| — | — | — | (B) All providers submitting specimens are informed of these criteria; |
| — | — | — | (C) All samples submitted are assessed for adequacy; |
| — | — | — | (D) Records of initial examinations and rescreening results are available; |
| — | — | — | (E) Rescreening of benign gynecological slides is: |
| — | — | — | (I) Performed by an individual who meets the personnel requirements for technical or general supervisor in cytology as defined under 42 CFR Part 493 Subpart M; |
| — | — | — | (II) Completed before reporting patient results on those selected cases; |
| — | — | — | (III) Performed and documented on: |
| — | — | — | · No less than ten percent of the benign gynecological slides; and |
| — | — | — | · Includes cases selected at random from the total caseload and from patients or groups of patients that are identified as having a high probability of developing cervical cancer, based on available patient information; |

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

- | | | | |
|---|---|---|---|
| — | — | — | (F) The technical supervisor: |
| | | | (I) Confirms all gynecological smears interpreted to be showing reactive or reparative changes, atypical squamous or glandular cells of undetermined significance, or to be in the premalignant (dysplasia, cervical intraepithelial neoplasia or all squamous intraepithelial neoplasia lesions including human papillomavirus-associated changes) or malignant category; |
| — | — | — | (II) Reviews all nongynecological cytological preparations; and |
| — | — | — | (III) Establishes, documents and reassesses, at least every six months, the workload limits for each cytotechnologist; |
| — | — | — | (G) All abnormal cytology reports are correlated with prior cytology reports and with histopathology reports if available, and the causes of any discrepancies are determined; |
| — | — | — | (H) Review of all normal or negative gynecological specimens received within the previous five years, if available in the laboratory system, or records of previous reviews, for each patient with a current high grade intraepithelial lesion or moderate dysplasia or CIN-2 or above; |
| — | — | — | (I) Notification of the patient's physician if significant discrepancies are found that would affect patient care and issuance of an amended report; |
| — | — | — | (J) An annual statistical evaluation of the number of cytology cases examined, number of specimens processed by specimen type, volume of patient cases reported by diagnosis, number of cases where cytology and histology are discrepant, number of cases where histology results were unavailable for comparison, and number of cases where rescreen of negative slides resulted in reclassification as abnormal; and |
| — | — | — | (K) Evaluation and documentation of the performance of each individual examining slides against the medical test site's overall statistical values, with documentation of any discrepancies, including reasons for the deviation and corrective action, if appropriate. |
| | | | (i) Immunohematology/ Transfusion Services: |
| — | — | — | (i) Perform ABO grouping, Rh (D) typing, antibody detection and identification, and compatibility testing as described by the Food and Drug Administration (FDA) under 21 CFR Part 606, and must also comply with 21 CFR Part 640. |
| | | | (A) Perform ABO grouping: |
| — | — | — | (I) By concurrently testing unknown red cells with FDA approved a anti-A and anti-B grouping sera; |
| — | — | — | (II) Confirm ABO grouping of unknown serum with known A1 and B red cells; |
| — | — | — | (B) Perform Rh (D) typing by testing unknown red cells with anti-D (anti-Rh) blood grouping serum; and |

Facility/MTS # _____ Date _____

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

— — — (C) Perform quality control procedures for immunohematology as described in Table 090-12.

Table 090-12 Quality Control Procedures - Immunohematology

Reagent	Control Material	Frequency
ABO antisera	• Positive control	• Each day of use
Rh antisera	• Positive and negative controls • Patient control to detect false positive Rh test results	• Each day of use • When required by the manufacturer
Other antisera	• Positive and negative controls	• Each day of use
ABO reagent red cells	• Positive control	• Each day of use
Antibody screening cells	• Positive control using at least one known antibody	• Each day of use

Y N NA WAC 246-338-090 (9) QUALITY CONTROL

— — — (ii) Blood and Blood Products:

— — — (A) Collecting, processing, and distributing:

— — — (I) Must comply with FDA requirements listed under 21 CFR Parts 606, 610.53, and 640; and

— — — (II) Must establish, document, and follow policies to ensure positive identification of a blood or blood product recipient.

— — — (B) Labeling and dating must comply with FDA requirements listed under 21 CFR 606, subpart G, and 610.53.

— — — (C) Storing:

— — — (I) There must be an adequate temperature alarm system that is regularly inspected.

— — — (II) The system must have an audible alarm system that monitors proper blood and blood product storage temperature over a twenty-four hour period.

— — — (III) High and low temperature checks of the alarm system must be documented.

— — — (D) Collection of heterologous or autologous blood products on-site:

— — — (I) Must register with the FDA; and

— — — (II) Have a current copy of the form FDA 2830 "Blood Establishment Registration and Product Listing".

— — — (iii) Must have an agreement approved by the director for procurement, transfer, and availability to receive products from outside entities.

— — — (iv) Promptly investigate transfusion reactions according to established procedures, and take any necessary remedial action.

Facility/MTS # _____ Date _____

Y N NA

WAC 246-338-090 (9) QUALITY CONTROL**(j) Histocompatibility:**

- | | | | |
|---|---|---|--|
| — | — | — | (i) Use applicable quality control standards for immunohematology, transfusion services, and diagnostic immunology as described in this chapter; and |
| — | — | — | (ii) Meet the standards for histocompatibility as listed in 42 CFR Part 493.1265, Condition: Histocompatibility, available from the department upon request. |

(k) Cytogenetics:

- | | | | |
|---|---|---|---|
| — | — | — | (i) Document: |
| — | — | — | (A) Number of metaphase chromosome spreads and cells counted and karyotyped; |
| — | — | — | (B) Number of chromosomes counted for each metaphase spread; |
| — | — | — | (C) Media used; |
| — | — | — | (D) Quality of banding; and |
| — | — | — | (E) Sufficient resolution to support the reported results; |
| — | — | — | (ii) Assure an adequate number of karyotypes are prepared for each patient according to the indication given for performing cytogenetics study; |
| — | — | — | (iii) Use an adequate patient identification system for: |
| — | — | — | (A) Patient specimens; |
| — | — | — | (B) Photographs, photographic negatives, or computer stored images of metaphase spreads and karyotypes; |
| — | — | — | (C) Slides; and |
| — | — | — | (D) Records; and |
| — | — | — | (iv) Perform confirmatory testing on all atypical results when performing determination of sex by X and Y chromatin counts. |

(l) Radiobioassay and Radioimmunoassay:

- | | | | |
|---|---|---|---|
| — | — | — | (i) Check the counting equipment for stability each day of use with radioactive standards or reference sources; and |
| — | — | — | (ii) Meet Washington State radiation standards described under chapter 70.98 RCW and chapters 246-220, 246-221, 246-222, 246-232, 246-233, 246-235, 246-239, 246-247, 246-249, and 246-254 WAC. |