



**WASHINGTON STATE WIC NUTRITION PROGRAM
HEMOCUE, INC. ORDER FORM**



WIC Agency Name: <hr/> Date Order Form Submitted: <hr/> Staff Completing Form: <hr/> Staff Phone Number: <hr/> Staff E-Mail Address: <hr/>	Clinic LASL ID Number: WAWIC - _____ Ship To: <hr/> <hr/> <hr/>
--	---

HemoCue Item #	Item	Number of Boxes to Order
121721	201+ HemoCue Analyzer – 1 per box	
111716	201+ Cuvettes - 200 per box (4 vials of 50)	
111715	201 + Cuvettes – 100 per box (4 small boxes of 25) / (Individually Wrapped)	
SLN200	Surgilance Lancet 1.8mm (100/Box)	
5252-58-104	Alcohol Prep Pads – 200 per box / (2” x 2” non-sterile, 2 ply, 70% isopropyl alcohol)	
1770-22082000	Gauze Pads – 200 per box / (2” x 2” 8 ply non-sterile)	
16-4823	Adhesive Bandages - 100 per box / (Large size: 3/4” x 3” sterile)	
668783	Adhesive Bandages - 100 per box / (Medium size: 5/8” x 2 1/4” sterile)	
1298033	Adhesive Bandages - 100 per box / (Small size: 3/8” x 1 1/2” sterile)	
139001	HemoCue Cleaners – 5 per pack/box	
564760	Sharps Containers – 1 per box / (Quart size)	
189011	Hard Shell Carry Case for 201+ HemoCue Analyzer – 1 per box	
CC ASSIGN	HB 201+ Control Cuvette	
22001002	Liquid Controls – HemoTrol – Low 1 ml vial (07.8 – 08.2 g/L) - 2 vials per box	If this is a first time ordering controls, prior to ordering, please contact Dick Court at 360-236-3694 to discuss your specific needs. The 201+ HemoCue Analyzer does its own self-calibration check and does not require controls. If your clinic is already on a delivery schedule for liquid controls and you need to make a change, please contact Patty Mattinson.
22002002	Liquid Controls – HemoTrol – Normal 1 ml vial (11.8 – 12.2 g/L) - 2 vials per box	
22003002	Liquid Controls – HemoTrol – High 1 ml vial (15.8 – 16.2 g/L) - 2 vials per box	

For any questions and to submit your HemoCue Order Form - Contact: Patty Mattinson - (360) 236-3606 Washington State WIC Nutrition Program P.O. Box 47886 Olympia, Washington 98504-7886 e-mail: patty.mattinson@doh.wa.gov / FAX: (360) 236-2320	<i>For use by HemoCue Only</i>
--	---------------------------------------

